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
API

Compact High Contrast Ratio Laser Attenuator
API

User Manual
Target version – v3.0.X

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- None

4.1.3. Connect()

Connects device via RS232 interface. USB adapter acts just as electrical adapter. If the device is already connected, this call returns false and keeps the old device connected.

Type	Name	Description
in <code>const std::string&</code>	<code>device_port_name</code>	Serial port. (e.g. "COM1")
return <code>bool</code>		Result of connection.

<Possible Error Code>

- None

4.1.4. Disconnect()

Disconnects device.

<Possible Error Code>

- None

4.1.5. GetCurrentStatus()

Queries device status.

Type	Name	Description
in <code>const int</code>	<code>serial_timeout</code>	Serial communication timeout in milliseconds (default=3500).

<Possible Error Code>

- AttenuatorStatus::OK
- AttenuatorStatus::SERIAL_XX
- AttenuatorStatus::MOTOR_XX

4.1.6. GetMotorInfo()

Queries motor information.

Type	Name	Description
out MotorInfo*	out_info	Return value.
in const int	serial_timeout	Serial communication timeout in milliseconds (default=7000).

<Possible Error Code>

- AttenuatorStatus::OK
- AttenuatorStatus::SERIAL_XX
- AttenuatorStatus::MOTOR_XX

4.1.7. MechanicalCalibration()

Start mechanical calibration. It should be called once after power is applied to the product.

Type	Name	Description
in const int	serial_timeout	Serial communication timeout in milliseconds (default=3500).

<Possible Error Code>

- AttenuatorStatus::OK
- AttenuatorStatus::SERIAL_XX
- AttenuatorStatus::MOTOR_XX

4.1.8. SetTransmittance()

Sets target transmittance value to run motor at the corresponding angle. The offset degree cannot be less than -360, and the sum of the offset degree and range degree cannot be more than 360.

Type	Name	Description
in const float	transmittance	Target transmittance value, range must be [0.0, 1.0].
in const double	offset_degree	Motor angle in degree when the transmittance is 0.

in	const double	range_degree	Angle in degree converted to transmittance range of 0-100% (default=45.0). That is, when the angle of the motor is offset_degree + range_degree, the transmittance will be maximum.
in	const int	serial_timeout	Serial communication timeout in milliseconds (default=3500).

<Possible Error Code>

- AttenuatorStatus::OK
- AttenuatorStatus::PARAMETER_OUT_OF_RANGE
- AttenuatorStatus::SERIAL_XX
- AttenuatorStatus::MOTOR_XX

4.1.9. GetTransmittance()

Queries current transmittance value. The offset degree cannot be less than -360, and the sum of the offset degree and range cannot be more than 360.

	Type	Name	Description
out	float*	out_transmittance	Return value, range might be [0.0, 1.0].
in	const double	offset_degree	Motor angle in degree when the transmittance is 0.
in	const double	range_degree	Angle in degree converted to transmittance range of 0-100% (default=45.0). That is, when the angle of the motor is offset_degree + range_degree, the transmittance will be maximum.
in	const int	serial_timeout	Serial communication timeout in milliseconds (default=3500).

<Possible Error Code>

- AttenuatorStatus::OK
- AttenuatorStatus::SERIAL_XX
- AttenuatorStatus::MOTOR_XX

4.2.0. SetMotorSteps()

Sets target steps to run motor at the corresponding angle.

Type	Name	Description
------	------	-------------

in	const int	steps	Target steps, range must be [-143360, 143360].
in	const int	serial_timeout	Serial communication timeout in milliseconds (default=3500).

<Possible Error Code>

- AttenuatorStatus::OK
- AttenuatorStatus::PARAMETER_OUT_OF_RANGE
- AttenuatorStatus::SERIAL_XX
- AttenuatorStatus::MOTOR_XX

4.2.1. GetMotorSteps()

Queries current steps.

Type	Name	Description
out int*	out_steps	Return value, range might be [-143360, 143360].
in const int	serial_timeout	Serial communication timeout in milliseconds (default=3500).

<Possible Error Code>

- AttenuatorStatus::OK
- AttenuatorStatus::SERIAL_XX
- AttenuatorStatus::MOTOR_XX

4.2.2. SetMotorDegree()

Sets target degree to run motor at the corresponding angle.

Type	Name	Description
in const float	degree	Target degree, range must be (-360.0, 360.0).
in const int	serial_timeout	Serial communication timeout in milliseconds (default=3500).

<Possible Error Code>

- AttenuatorStatus::OK
- AttenuatorStatus::PARAMETER_OUT_OF_RANGE
- AttenuatorStatus::SERIAL_XX
- AttenuatorStatus::MOTOR_XX

4.2.3. GetMotorDegree()

Queries current rotary angle of the motor.

Type	Name	Description
out float*	out_degree	Return value, range might be [-143360, 143360].
in const int	serial_timeout	Serial communication timeout in milliseconds (default=3500).

<Possible Error Code>

- AttenuatorStatus::OK
- AttenuatorStatus::SERIAL_XX
- AttenuatorStatus::MOTOR_XX

4.2.4. StartOptimization()

Optimizes the operating frequencies for backward and forward movement by searching frequencies and monitoring mechanical performances. Allowing 20 minutes cool down period is highly recommended after optimization. This operation may take several minutes, during that time the device and the associated communication bus are blocked. In this case, AttenuatorStatus::MOTOR_BUSY is thrown when another method is called, except Stop().

Type	Name	Description
in const int	serial_timeout	Serial communication timeout in milliseconds (default=3500).

<Possible Error Code>

- AttenuatorStatus::SERIAL_XX
- AttenuatorStatus::MOTOR_XX

4.2.5. StartCleaningMechanics()

Starts cleaning cycle, during which the device will move backwards and forwards over its full range of travel for a few minutes, in order to remove any dust or debris from the bearing rails and motor contacts. Allowing 20 minutes cool down period is highly recommended after optimization. This operation may take several minutes, during that time the device and the associated communication bus are blocked. In this case, AttenuatorStatus::MOTOR_BUSY is thrown when another method is called, except Stop().

Type	Name	Description
in const int	serial_timeout	Serial communication timeout in milliseconds (default=3500).

<Possible Error Code>

- AttenuatorStatus::SERIAL_XX
- AttenuatorStatus::MOTOR_XX

4.2.6. Stop()

Aborts the optimization or cleaning process initiated by StartOptimization() or StartCleaningMechanics().

Type	Name	Description
in const int	serial_timeout	Serial communication timeout in milliseconds (default=3500).

<Possible Error Code>

- AttenuatorStatus::OK
- AttenuatorStatus::SERIAL_XX
- AttenuatorStatus::MOTOR_XX

4.2. HILLAB::AttenuatorStatus Enum

Value	Description
PARAMETER_OUT_OF_RANGE	User input error
SERIAL_CLOSED	Unable to write bytes to serial port
SERIAL_TIMEOUT	Unable to receive bytes from serial port
SERIAL_UNDEFINED_ERROR	Unexpected serial error
OK	No problem
MOTOR_COMMUNICATION_TIMEOUT	Device internal error
MOTOR_MECHANICAL_TIMEOUT	Device internal error
MOTOR_COMMAND_ERROR	Device internal error
MOTOR_VALUE_OUT_OF_RANGE	Device internal error
MOTOR_MODULE_ISOLATED	Device internal error
MOTOR_MODULE_OUT_OF_ISOLATION	Device internal error
MOTOR_INITIALIZING_ERROR	Device internal error
MOTOR_THERMAL_ERROR	Device internal error
MOTOR_BUSY	Device internal error
MOTOR_SENSOR_ERROR	Device internal error
MOTOR_MOTOR_ERROR	Device internal error
MOTOR_OUT_OF_RANGE	Device internal error
MOTOR_OVER_CURRENT_ERROR	Device internal error

4.3. HILLAB::MotorInfo Struct

Type	Name	Description
std::string	serial_number	Motor serial number
std::string	year_of_manufacturing	Motor year of manufacturing
std::string	travel_degree	Motor travel degree
std::string	pulse_per_degree	Motor pulse per degree
std::string	firmware_release	Motor firmware release information
std::string	hardware_release	Motor hardware release information
std::string	mechanical_offset	Motor factory offset

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