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

The screenshot displays the HILLAB Settings application window. The interface includes a menu bar (File, View, Devices, Help) and a toolbar with 'Run Stream' and 'Settings...' buttons. The 'Settings' window is open, showing a sidebar with categories: Device, Noise Reduction, Profiling, Visualization, Server, and General. The 'Device' tab is selected, displaying various configuration options. Red circles with numbers 1 through 11 highlight specific controls:


- 1: HILLAB logo
- 2: Device Alias input field
- 3: Acquisition Trigger section, including 'Auto' and 'External Trigger' checkboxes
- 4: FPS slider control
- 5: Trigger Activation dropdown menu (set to 'Rising Edge')
- 6: Exposure slider control
- 7: Auto Exposure Once button
- 8: Gain slider control
- 9: Auto Gain Once button
- 10: Strobe Output toggle switch
- 11: 2.5V Voltage Output toggle switch

The right side of the image shows a zoomed-in view of the settings controls, highlighting the Exposure, Auto Exposure Once, Gain, Auto Gain Once, Strobe Output, and 2.5V Voltage Output sections.

Index	Key	Type	Method	Value
12	profile.lowpassFilterEnabled	bool	Get: 값 가져오기 Set: 값 변경	true = 활성화 / false = 비활성화
13	profile.backgroundRemovalEnabled	bool		true = 활성화 / false = 비활성화
14	profile.backgroundRemovalOffset	float		범위: 0.001~1.000 / 단위: 0.001
15	profile.medianFilterEnabled	bool		true = 활성화 / false = 비활성화
16	profile.medianFilterKernelSize	int32		범위: 3~21 / 조건: 홀수
17	profile.thresholdEnabled	bool		true = 활성화 / false = 비활성화
18	profile.thresholdGaussianSigma	float		범위: 0.0~3.0 / 단위: 0.1
19	profile.thresholdFactor	float		범위: 0.000~1.000 / 단위: 0.001
20	profile.clipLevel	float		범위: 0.000~1.000 / 단위: 0.001
21	profile.averageCount	int32		범위: 1~10 / 단위: 1
22	profile.targetPointX	float		단위: 0.1
23	profile.targetPointY	float		단위: 0.1
24	profile.dataAnalysisPeriod	int32		범위: 3~21 / 조건: 홀수
25	profile.pixelToMetric	bool		true = 측정 단위 → 마이크로미터 false = 측정 단위 → 픽셀
26	profile.pixelToMetricScale	float		범위: 0.000000~100.000000 / 단위: 0.000001

Index	Key	Type	Method	Value
27	image.colormap	int32	Get: 값 가져오기 Set: 값 변경	0 = None 1 = Jet 2 = Turbo 3 = Rainbow-Short 4 = Rainbow-Long 5 = Thermal 6 = Cool-Warm 7 = Viridis
28	image.highlightIntensityLimits	bool		true = 시각화 활성화 false = 시각화 비활성화
29	image.drawCentroid	bool		
30	image.drawTargetPoint	bool		
31	image.drawDiameter	bool		0 = D4σ 1 = FWHM 2 = 1/e <sup>2</sup>
32	image.diameterType	int		
33	image.drawIntensity	bool		true = 시각화 활성화
34	image.drawUnit	bool		false = 시각화 비활성화

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